

# 48 hours of hydration skinsolver<sup>®</sup> protocol

## step 1 (1 minute)

Begin the Dermalogica Double Cleanse by applying 2 pumps of PreCleanse with fingers, cotton pad or fan masque brush. Swipe over skin to remove make-up, oil and debris.

## step 2 (3 minutes)

Mix 1 pump Special Cleansing Gel or skin specific cleanser with water in a cup and apply with fan masque brush. Lightly manipulate, then remove with cool damp disposable wipes or cool steam towel.

**Pro Product Boost:** Add 1 pump of Exfoliant Accelerator 35 to the cleanser.

**Modality Option:** Use an ultrasonic brush or blade on exfoliation mode for a deep cleanse.

## step 3 (3 minutes)

Dispense a few drops of Skin Hydrating Booster and press onto the skin or add into the masque. Then apply a thin layer of Skin Hydrating Masque with a fan masque brush. Remove with damp disposable wipes.

**Modality Option:** Use an ultrasonic blade on penetration mode to optimise hydration.

## step 4 (2 minutes)

Spritz the skin with Multi-Active Toner and apply Skin Smoothing Cream.

## step 5 (1 minute)

Protect skin by layering Solar Defense Booster SPF50 or prescribed DaylightDefense.

